

Colorwatcher, Colorimetric Array Reader



Colorwatcher is a compact and portable instrument designed for imaging colorimetric arrays (for example, DNA, protein biochips). The instrument by standard operates in brightfield and darkfield illumination modes with bottom and top illuminations respectively. The instrument is able to work with 96-well plate, a 12 x 8 well strip, as well as microscopic glasses (placed in a special adapter). The device is equipped with a highly sensitive CMOS camera. The instrument has an automatic mechanism of 96 well plate ejection outside the body.

The instrument can operate with a 20V portable lithium-ion battery, which makes the device usable in the field studies or point-of-care.

The compact design allows the device to be carried in the hand luggage of the aircraft.

Software Features:

- Automatic array finder via machine learning and image recognition
- Grid lay-outing (Manual and Automatic)
- Analysis of images by the average/median intensity of the spots
- Qualitative/quantitative analysis of the arrays
- Creating Qualitative/quantitative analysis assays
- Quantitative assays with 4/5 parameter logistics functions, etc.
- Setting multi level interpretation thresholds for different type of samples in the same well (e.g. tolerance to egg and lettuce)
- Report exports to PDF, CSV, EXCEL
- Control of camera exposure, gain, XYZ kinematics.

Colorwatcher, Colorimetric Array Reader

Specifications	
Illumination	Brightfield and darkfieled (bottom and top)
Vessels	96 well plate/12 x 8 well strip/4 microscope slides
Light source	White LED
Lifetime of the light source	ce <10 000 hours
Data interface for unit co	ntrols USB 2.0
Camera	3 MP, CMOS
Resolution	5-7 μm per pixel, per user request
Data interface for camera	control USB 3.0
Magnification	>10x
Image formats	png/bmp/tiff/jpg, 8/16 bit
Focus	Manual, adjustable via PC
Software	Included
PC requirements	CPU: Intel i7 (with AVX2 support) RAM: 8 GB Video card: Nvidia GTX 1050 Ti 4GB, or better (Capability only with Nvidia cards) SSD: 256 GB OS: Windows 10 (64 bit)
Size (W×D×H)	330 × 345 × 150 mm
Weight, w/o power suppl	y not more than 6 kg

Power supply

Input AC 100–240 V 50/60 Hz, Output DC 20-24 V, 2.5A

Software interface, picture of a well of 96 well plate (brightfield)

it controls Image analyser	Data analyser																											
een Red Load																												
IR BIR CIR DIR	E1 R F1 R	G1 R H1 R	A2 R B	2 R C2	R D2	R E2	R F2	R G2 R	H2 R	A3 R	B3 R C3	R D3 F	E3 R F3	R G3 R	H3 R	A4 R E	4 R 0	1 R D4	R E4 F	F4 R	G4 R	H4 R	A5 R	85 R C5 R	D5 R E	5 R F5 R G5 R H5	R A6 R	86 R C
Manual Auto				-		-			_						Input D	sta		View Results										
Manual controls:					-											2	3	4	5	6	7	8	9	10		Choose an assay		
C Enable Grid										1				1	1	2	3	4	3	U		0	9	10		· · · · ·		ASSAY
15 🔹 LEN			/											-	54555	19209	6686	18715	6631	6692	6684	6661	6729	6624		Choose a Template or S	iave as	
9 🖨 SHIFT											1			2												1 What to show in a	a cel:	
180 SHIFT_X1															19299	6688	6616	7127	6820	7082	6597	6830	6467	6599		A ////// Sample Name	~	
118 SHIFT_Y1														3												RFU	~	
180 SHIFT_X2				•	• •	0 0	0 0	0 0	0 0	0		/			7280	6724	6729	23823	6809	7311	24899	6560	6572	6407		Cell Name Sample Name	0	Type
118 SHIFT_Y2				C	0	0 0	0 0	o c	0 0	0		1		4												Man cha	nal .	
-60 SHIFT_X				C	0 0	0 •	0 0	• 0	0 0	0					34170	6856	6605	6746	6653	6572	51969	6533	6402	6539		Calculate Highlight	~	
-47 SHIFT_Y					• •		0 0	• •	0 0	0				5														
121 SHIFT_X_B				C	0	0 0	0 0	0 0	• •	0					7788	16945	7001	6613	6613	6551	6447	43858	6255	6192				
56 SHIFT_Y_B				C	0 0	0 0	0 0	o c	0 0	0				6						1000012			2000017	10000				
0.00 C Rotate						0 0	•	o c	0 0	0				-	19437	18608	6793	27209	6588	6564	6372	6403	6274	6555				
				C	0	• •	0	ò (0 0	•			/	7	55017	29966	6640	10692	25512	7244	6463	6346	6407	8359				
Calculate				c	0 0	0 0	0	0 0	0 0			1		8	33917	29900	0048	10092	55512	7344	0403	0340	0407	8339				
Calculate All Clone to next		1			0		0	0 0	0 0	•		/			16516	6809	46475	6643	6825	36223	6868	6475	10672	23945				
Clone to all				-						-				9														
Clone to all remaining			2.												17078	16136	6527	7637	6623	6640	6631	6672	29237	33811				
Load image											/			1														
Save image														0	55471	6766	6801	7163	7239	6772	6671	6704	6294	55865				
Save Coordinates										/																		
Save All Coordinates																												
	Screenshot	Туре		Gree	n			-	Red				Measure	i L														
		• 96 Well		1.1.1.	1.0	12 12	1-1				_	_	Go To	i —														
		Oslide		2 3	4	5 6	7	8 9	10 1	1 12	✓ Ched	k All	Stop															
		Uside	^ 🖌		∕∠∠		∕∠∠		106		\checkmark		otop	1														
		CloneTo	B 🗸		70		70		100		\checkmark		Reference															
			· 7		70		70		100	77			Home	1														
					ลักั		žĂ	Ĩ		лй.																		
			E										Eject															
								<u>N</u>		39																		
	State				<u>í</u>		<u>1</u>																					
	Additional		<u>د</u>		2⊻		⊴⊻																					
	Com panel		н 🔽		∕∕∕		∕∕∕		10																			
	Preset:											- 11																